

Mass vs. Diameter - A new approach for analysing particle size and structure

Franz Friebe, Jörg Wieder, Kingsley Reavell
femtoG AG, Winterthurerstrasse 292, 8057 Zurich, franz.friebe@femtoG.com

What is the meaning of the diameter for non-spherical particles?

Nanoparticle – named after a length, but traded in €/kg, grams, big bags, ...




How many nanoparticles are in one gram of powder?

Analysis of particle mass and diameter

→ structure of engineered nanomaterials.

8-12 min per scan; 10-50 mio particles counted



poster download





reach out to me

Carbon Black

Increase in aggregation level at constant PP-diameter of 27nm
→ Increase in mass and size
→ Decrease in density and aggregates/gram

$$V_{\text{Pore}} \approx \text{oil number} \approx \text{density}^{-1}$$



Inorganic Pigments



What is the true particle size?

Results from most sizing methods are affected by different refractive indices and structures

Median mass	Particles per gram	«Fractal»-index
Cr ₂ O ₃ 350nm sphere @ 5.2 g/cm ³ m _{sphere} : 117fg → isolated spherical particles	4.8 · 10 ¹² /g	2.85
Fe ₃ O ₄ - 300nm sphere @ 4.6 g/cm ³ m _{sphere} : 65 fg → isolated spherical particles	1.1 · 10 ¹³ /g	2.80
Fe ₂ O ₃ - 170nm sphere @ 5.0 g/cm ³ m _{sphere} : 13 fg → low aggregation level	9.2 · 10 ¹² /g	2.88
Fe ₂ O ₃ - 90nm sphere @ 5.0 g/cm ³ m _{sphere} : 1.9 fg → fractal-like aggregate	1.8 · 10 ¹³ /g	2.72
FeOOH rods: 100 x 400nm @ 4.0 g/cm ³ m _{rod} : 16fg	9.5 · 10 ¹² /g	2.70

Two-component system,
→ Optical sizing methods fail

Coated and stabilized pigments



Median mass	Particles per gram	Mass fraction		«Fractal»-index
		measured	spec-sheet	
AlPO ₄ isolated non-aggregated particles				
3.8 fg	1.1 · 10 ¹³ /g	10%	26%	2.98
FeOOH rods: 100 x 700nm @ 4.0 g/cm ³ (+coating) m _{rod} : 29fg (+7fg)				
115 fg	4.7 · 10 ¹² /g	90%	(74%)	2.56

Only 61% of the AlPO₄ is used for coating.
→ Reduced color strength, loss of material, ...

1. Sampling

a) from reactor:

- dilution
- cooling
- drying

b) powder

- dry-dispersion
- solvent-spray dispersion


2. Deagglomeration

- ultrasound
- venturi-nozzles
- adjustable deagglomeration intensity
- shear: 100-500 N/m²
- force: 5 · 10⁻¹⁰ N / particle

probing particle stability

3. Centrifugal particle mass analyzer

- selection by mass to charge ratio



4. Particle sizing

- mobility to size ratio

